


The magnetic stirring action of «alluro» (MRK 2): Power is transmitted from the stirrer motor to the drive collar by means of a flexible coupling. The drive collar is fitted with a magnetic ring, which uses the mutual attraction of the magnetic fields to couple with the inner magnet. This means that the stirrer shaft is driven without any contact.



alluro magnetic stirrer head



Torque	20 Ncm	40 Ncm	60 Ncm	90 Ncm	20 Ncm	40 Ncm	60 Ncm	90 Ncm
Product no.	2/20	2/40	2/60	2/90	2/20 HC22	2/40 HC22	2/60 HC22	2/90 HC22
Material	WNr. 1.4435	WNr. 1.4435	WNr. 1.4435	WNr. 1.4435	HC22	HC22	HC22	HC22
WNr./AISI	AISI 316 L	AISI 316 L	AISI 316 L	AISI 316 L	-	-	-	-
T °C	240	240	240	240	240	240	240	240
Max. vol. ml	2000	4000	6000	10000	2000	4000	6000	10000
A mm	44	56	78	100	44	56	78	100
B mm	156	168	190	212	156	168	190	212
C mm	113	125	147	169	113	125	147	169

Connection and seal

- Ground joint NS 45/40
- Double seal to glass reactor, with two consecutive Viton O rings. In the Hastelloy or titanium versions, the grooves increase in size, and the quality of the O rings changes to EPDM or Kalrez.

Material

- All parts coming into contact with the medium are made of corrosion-resistant material DIN 1.4435. We offer the same line in Hastelloy C22 for aggressive media. Titanium, tantalum and types of Inconel can also be used.

Drive / torque and bearing

- The best solution is to connect the MRK to a factory-fitted drive motor via a flexible stirrer component. The standard drive shaft of the MRK is 7 mm in diameter. A square 6 mm connection is also available.
- The torque of the magnetic coupling can be increased from 20-90 Ncm.
- Corrosion-resistant ball bearings are used in the driven shaft for a maximum speed of 3000 rpm. Hastelloy or titanium versions use friction bearings for a maximum speed of 1600 rpm.

delivery address

correspondence address

phone

fax

internet

e-mail

premex reactor ag

industriestrasse 11

p.o. box 444

2543 lengnau/switzerland


+41 (0)32 653 01 52

+41 (0)32 652 11 80

www.premex-reactorag.ch

office@premex-reactorag.ch

omni
mrk



premex
reactor ag 